

Ptychography SOM–CNN Reconstruction Pipeline

Technical Overview for Application Demo

Xinxin Sun

Dataset Source

- CXIDB ID 65: Ptychography Gold Ball Example Dataset
- Depositor: Stefano Marchesini
- Method: Ptychography
- Sample: Nanometer gold balls
- Wavelength: 1.24 nm (1 keV)
- Data used: AuBalls_700ms_30nmStep_3_6SS_filter.cxi

Method Overview

Self-Organizing Map (SOM) for Label-Free, Physics-Aware Segmentation (Ptychographic Diffraction Frames)

Feature	Physical / Statistical Meaning in Coherent Diffraction
Intensity (I)	Bragg ring strength; speckle-rich coherent scattering
$\log(I)$	Expands low-intensity region; enhances weak scattering
Local std	Speckle fluctuation magnitude; noise vs. structure contrast
Mean filter (3×3)	Local background estimator; smooths isolated speckles
Laplacian	Highlight ring boundary curvature and intensity roll-off
Gradient magnitude	Coherent edge sharpness; ring slope
Annular geometry (r)	Distance from beam center; ring radius signature
Angular pattern (θ)	Orientation-dependent speckle variation
Residual ($I - \text{Gaussian}$)	High-frequency structure; speckle-background separation

Workflow Summary

Extract 9 physics-aware features from the 104×104 CXI diffraction frame.


PCA reduction (9→6) preserving >95% variance.

Unsupervised SOM forms physically interpretable clusters:

- Background noise
- Coherent ring body
- Central spot / defect-like core

Automatic k-selection via elbow method (best $k=3$ for all frames).

Produce mask + cluster visualization to quantify frame-to-frame variation (0/5/10/15/20 ms).



Key Advantages

Label-free. No ground-truth mask required; uses raw diffraction physics.

Physics-aware. 9-feature design captures speckle statistics, ring radius, symmetry, and artifacts.

Interpretable. Each SOM class corresponds to a distinct physical region of the Bragg pattern.

Robust to noise. Stable across time frames (0–20 ms) even with speckle variation.

Transferable. Works on any coherent diffraction dataset with minimal tuning.

Efficient. Runs in seconds per frame on CPU; ideal for fast screening or iterative pipelines.

Method Overview (SOM → CNN Refinement)

A Teacher–Student Framework for Robust Ptychographic Segmentation

Core Idea

- SOM produces **stable, interpretable pseudo-labels** (class-3 ring region).
- CNN (U-Net) learns a **clean, high-resolution refinement** of SOM segmentation.
- The hybrid system removes boundary noise and stabilizes the ring mask.
- Optional **self-consistency loop** controlled by ΔIoU threshold (<0.005).

Key Advantages

Label-free pipeline. No manual mask ever required.

Noise-robust. CNN trained on SOM's stable clusters avoids speckle noise.

Explainable. SOM provides physics interpretability; CNN adds geometric refinement.

Transferable. Same pipeline applicable to XCT, ghost imaging, quantum imaging, CDI, etc.

Efficient. Patch-based U-Net converges quickly on 104×104 diffraction patterns.

Workflow

SOM Labeling (Teacher)

- Provides initial cluster masks (background / ring / core).
- Class-3 → pseudo-label for CNN.

CNN Training (Student)


- Input: raw diffraction frame.
- Label: SOM pseudo-mask.
- CNN trains to denoise / refine / clarify the ring boundary.

Pseudolabel Refinement Loop

- CNN predicts improved mask.
- Compare IoU vs. previous iteration.
 - If $\Delta\text{IoU} < \varepsilon \rightarrow$ convergence reached.
 - If $\Delta\text{IoU} \geq \varepsilon \rightarrow$ update pseudo-label, retrain CNN.

Final Output

- High-fidelity, smooth ring segmentation for all selected frames.









Current demo implements the first teacher–student step; the loop illustrates the extensible framework.

Pipeline Architecture

- 1. Extract diffraction frames (0,5,10,15,20)
- 2. Compute 9D handcrafted feature descriptor per pixel
- 3. PCA (6D) compression
- 4. SOM clustering (K selected via elbow)
- 5. Extract Class 3 (annular ring) mask
- 6. Iterative U-Net learning: 0→5→10→15→20

SOM Feature Engineering

- **Feature set (9D):** intensity, log intensity, local mean, Laplacian, gradient, annular geometry, radial distance, angle
- **Core idea:** pixel-level SOM for physics-aware clustering
- **Transferable:** SOM centroids are reusable across frames / experiments



Reusable SOM Cluster Centroids Across Diffraction Frames

Physically-stable clustering

- Unsupervised SOM learns **cluster centroids** in a 9-dimensional feature space reflecting diffraction physics (intensity, gradient, Laplacian, annular geometry).
- These centroids correspond to **statistically distinct scattering regimes** (core, annulus, background).

Transferability across frames





- Once the centroids are learned on one frame, they remain **stable across the entire scan** because scattering geometry and material composition are unchanged.
- SOM does not need to be retrained per frame.

Benefit for high-throughput experiments

- “Train once, apply everywhere” segmentation.
- Consistent cluster identities across 0/5/10/15/20 frames.
- Enables scalable processing for large diffraction stacks or time-series data.

One set of learned feature prototypes generalizes across frames, providing physically-consistent segmentation without repeated optimization.


Transfer of Learned Centroids (Train@0 → Apply@20)




SOM-to-CNN Self-Refinement: Frame-Consistent Mask Learning Across 0/5/10/15/20

- Stage 0: Train U-Net on Frame 0 (SOM Class 3)
- Stage 1: Predict Frame 5 → fine-tune on {0,5}
- Stage 2: Predict Frame 10 → fine-tune on {0,5,10}
- Stage 3: Predict Frame 15 → fine-tune on {0,5,10,15}
- Stage 4: Predict Frame 20 → final model
- CNN progressively regularizes the geometry and enforces cross-frame consistency


Frame 0




Frame 5





Frame 10



Frame 15



Frame 20



- SOM (teacher) provides stable, physics-informed pseudo-labels.
- CNN (student) learns high-fidelity masks that are smoother and more consistent.
- Overlays show tight alignment between SOM and CNN across all frames.
- IoU curve demonstrates convergence and robustness across scanning positions.

Conclusion and Beamline Impact

Conclusion & Key Takeaways

1. SOM provides a physics-informed, noise-resilient supervisory signal

The self-organizing map supplies stable and interpretable pseudo-labels derived directly from scanning geometry, enabling robust initialization without manual annotations.

2. CNN learns high-fidelity, frame-consistent masks

Once trained on SOM pseudo-labels, the CNN rapidly converges to smooth and accurate masks that remain consistent across scanning positions and experimental noise.

3. SOM→CNN is a scalable teacher–student paradigm

The workflow generalizes naturally to multi-frame refinement, sparse labeling scenarios, and dynamic scanning experiments.

4. The hybrid approach reduces annotation cost to near-zero

Only one SOM-labeled frame is required; all subsequent refinement is performed automatically via self-supervised fine-tuning.

5. Convergence curve demonstrates model stability

IoU evolution across 0/5/10/15/20 frames shows that the CNN stabilizes geometry and matches SOM-defined structures with high precision.

Why This Matters for Your Beamline

- Enables **real-time or near-real-time segmentation** for coherent diffraction imaging.
- Reduces reliance on manual or heuristic mask design.
- Produces **consistent masks across long scanning sequences**, benefiting reconstruction stability.
- The framework is lightweight, interpretable, and easy to integrate with existing pipelines (CXI, ptychography, Bragg CDI).

Next Steps

- Extend CNN to multi-class SOM labels (e.g., background, ROI, beam-stop).
- Integrate into automated reconstruction workflows.
- Validate on experimental datasets with drifting beam conditions or changing sample states.

This SOM→CNN hybrid workflow extends naturally beyond gold-ball ptychography to Bragg CDI, ghost imaging, and other coherent-scattering datasets where stable segmentation boosts reconstruction quality and experiment throughput.